

PEI Transfection

Transfection:

1. Split 293T cells one day before transfection in DMEM/10% FBS medium:
 - a. **6 well dish:** 0.5×10^6 cells
 - b. **10cm dish:** 4.0×10^6 cells
 - c. **15cm dish:** 9.0×10^6 cells
2. Prior to transfection bring all reagents to room temperature.
3. In a sterile tube dilute total plasmid DNA (ug) in **serum-free** DMEM w/o phenol red (volume of media is 10% of final volume in culture vessel). Use transgene: viral packaging (psPAX2):viral envelope (pMD2G) constructs at 4:2:1 DNA ratio
 - a. **6 well dish:** 200ul + 3 ug of total DNA
 - b. **10cm dish:** 1mL + 7-8 ug of total DNA
 - c. **15cm dish:** 2mL + 11-12 ug of total DNA
4. Add PEI (1ug/uL) to the diluted DNA. Mix immediately by vortexing or pipeting. The volume of PEI used is based on a 3:1 ratio of PEI (ug):total DNA (ug).
 - a. **6 well dish:** 9ul of PEI(1ug/ul) = 9ug
 - a. **10cm dish:** 21ul of PEI (1ug/ul) = 21ug
 - b. **15cm dish:** 33ul of PEI(1ug/ul) = 33ug
5. Incubate 15 minutes at RT
6. Add DNA/PEI mixture to cells
7. Harvest transfected cells and/or viral supernatant at 48 hours post-transfection

Helper	6ug
AAV9	4ug
crispr	4ug

PEI: 42ul

Reagents:

PEI (1ug/ul) – PEI is Polyethylenimine 25kD linear from Polysciences (cat# 23966-2).

To make a stock solution:

- Dissolve PEI in endotoxin-free dH₂O that has been heated to ~80°C.
- Let cool to room temperature.
- Neutralize to pH 7.0, filter sterilize (0.22um), aliquot and store at -20°C; a working stock can be kept at 4°C.